



东莞中力仪器科技有限公司

Dongguan Zhongli Instrument Technology Co., Ltd

Tel : +86-769-33657986

Fax: +86-769-33657987

Website : [www.zltester.com](http://www.zltester.com)

E-mail: [Zhongli@zltester.com](mailto:Zhongli@zltester.com)

Add: No.62, Xinhe Management of Area, Wanjiang District, Dongguan, Guangdong, China

# Introduction

**Name: ZL-6014B Desktop Electric Oven (Blast Drying Oven)**



## Equipment summary:

This machine uses high-temperature to dry heat, it also can make micro-organisms oxidation, protein denaturation, to kill residual microorganisms in a certain period.

## Technical parameter:

Model	6014B- I	6014B- II	6014B-III	6014B- V	6014B-VI
Power	220V 50Hz				
temperature	50℃—200℃				



# 东莞中力仪器科技有限公司

Dongguan Zhongli Instrument Technology Co., Ltd

Tel : [+86-769-33657986](tel:+86-769-33657986)

Fax: [+86-769-33657987](tel:+86-769-33657987)

Website : [www.zltester.com](http://www.zltester.com)

E-mail: [Zhongli@zltester.com](mailto:Zhongli@zltester.com)

Add: No.62, Xinhe Management of Area, Wanjiang District, Dongguan, Guangdong, China

fluctuation	±1℃				
Following alarm	+5℃				
Consume power	500W	500W	800W	1500W	2000W
Inside size (L*W*H) mm	250×250×250	300×300×275	415×370×345	554×370×523	600×550×600
Outside size (L*W*H) mm	535×440×435	585×490×460	700×560×530	840×560×710	885×740×790

## Main function:

This Blast Drying Oven can be used in drying glass, the thermal denaturation of sample and thermal in-duration experiment, eliminate specimen's residual moisture, microbes in the vessels of the glass with dry heat sterilization, experiment of preheating, before heating food and chemical need to harden or softening experiment, ruled out moisture, extend the food in the harsh conditions of the-preservation time, dry and heat resistance experiment the building materials, electronic component-sand dry, baking, disinfection, sterilization and material degeneration test.

## Main feature:

1. The vertical drum wind drying oven is made of high qualified steel plate or mirror stainless steel, plastic spraying, studio shell double-layer tempered glass windows, beautiful model-ling novel.
2. ZL-6014B vertical electric drum wind microcomputer intelligent control, temperature drying oven overheating protection, time settings, temperature, small fluctuations linear
3. ZL-6014B series vertical drum wind drying oven instrument has tracking alarm function of mixing hot air circulation, the heat loss is small, the drum wind, heat can be independent operation.